

Data Sheet

Model: BCF12/25-SY

General Purpose Chamber Furnace



INTRODUCTION

Sylab Scientific Pvt. Ltd., under technical collaboration with Elite Thermal Systems Ltd, UK, manufactures high quality furnaces in India.

Elite Thermal Systems, based in UK, is a 30-year-old furnace manufacturing company with strong brand recognition across the globe.

Sylab Scientific delivers high-quality solutions by leveraging technical expertise and components sourced from Elite Thermal Systems Ltd, UK.

BCF12-SY General purpose chamber furnace is designed for General purpose laboratory applications and providing satisfactory performance for many firing applications and is available in various capacities ranging from 5 litres to 45 litres in standard versions. Other chamber capacities/volumes are offered in custom-built models.

The **BCF12/25-SY** chamber furnace has a 25-liter capacity.



BCF12/25-SY

SPECIFICATIONS

Maximum Temperature: 1200°C

Maximum Continuous Temperature: 1150°C

- | Chamber dimensions (mm): 250 x 250 x 400 (H x W x D)
- | Two-sided panel heating elements of ceramic fibre containing wire spirals freely radiating from sinusoidal grooves
- | Compact bench mounted design
- | Energy efficient, high quality, low thermal mass insulation
- | Vertically lifting door keeps the hot surface away from the user
- | Positive break door safety switch isolates heating elements from power supply when door is open
- | ELITE PID controller offering 8 programs with 8 segments each
- | Small ceramic chimney and Hard Ceramic hearth tile are fitted as standard
- | **External Dimensions (mm):** 780 x 635 x 695 (H x W x D) (Indicative)
- | **Net Wt.:** 67 kg
- | **Supply / Power:** 230V– 1 Phase – 6.0 kW

OPTIONS

- | Independent Over-temperature protection system
- | Multi segment, multi program storage Controller
- | Powered exhaust, Chimney with fan

SYLAB SCIENTIFIC PVT. LTD.

Formerly SYLAB PVT. LTD.

Ground Floor, B-50, Industrial Estate, Sanath Nagar, Hyderabad-500018. INDIA,
Phone: +91 40 67216376 E-mail: sales@sylabscientific.com Website: www.sylabscientific.com