

Data Sheet

Model: BCF13/5-SY

General Purpose Chamber Furnace



INTRODUCTION

Sylab Scientific Pvt. Ltd., under technical collaboration with Elite Thermal Systems Ltd, UK, manufactures high quality furnaces in India.

Elite Thermal Systems, based in UK, is a 30-year-old furnace manufacturing company with strong brand recognition across the globe.

Sylab Scientific delivers high-quality solutions by leveraging technical expertise and components sourced from Elite Thermal Systems Ltd, UK.

BCF13-SY is a fast-heating furnace for general purpose use, where clean operating conditions prevail. Good temperature uniformity is achieved by the use of open heating elements retained in low thermal mass chamber wall panels with capacity ranging from 5-litre to 45-litre in standard versions. Other chamber capacities/volumes are offered in custom-built models.

The **BCF13/5-SY** chamber furnace is a fast-heating furnace with 5 litre capacities.



BCF13/5-SY

SPECIFICATIONS

Maximum Temperature: 1300°C

Maximum Continuous Temperature: 1250°C

| **Chamber dimensions (mm):** 150 x 150 x 225 (H x W x D)

| **Designed for** light duty laboratory applications

| **Energy efficient,** high quality, low thermal mass insulation

| **Two-sided panel heating elements** of ceramic fibre containing wire spirals freely radiating from sinusoidal grooves

| **Insulated door opens upwards & outwards** keeping the hot face away from the operator when the door is opened

| **Vertically lifting door** keeps the hot surface away from the user

| **Positive break door safety switch** isolates heating elements from power supply when door is open

| **ELITE PID controller** offering 8 programs with 8 segments each

| **Small ceramic chimney** is fitted as standard

| **External Dimensions (mm):** 593 x 400 x 505 (H x W x D)
(Indicative)

| **Net Wt.:** 36 kg

| **Supply / Power:** 230V– 1 Phase – 2.0 kW

OPTIONS

| **Independent Over-temperature protection system**

| **Multi segment, multi program storage Controller**

| **Powered exhaust, Chimney with fan**

SYLAB SCIENTIFIC PVT. LTD.

Formerly SYLAB PVT. LTD.

Ground Floor, B-50, Industrial Estate, Sanath Nagar, Hyderabad-500018. INDIA,
Phone: +91 40 67216376 E-mail: sales@sylabscientific.com Website: www.sylabscientific.com