

Data Sheet

Model: TMH15(Z3)/50/610-SY

1500°C Maximum Temperature,
Three Zone Horizontal Tube Furnace



INTRODUCTION

Sylab Scientific Pvt. Ltd., under technical collaboration with Elite Thermal Systems Ltd, UK, manufactures high quality furnaces in India.

Elite Thermal Systems, based in UK, is a 30-year-old furnace manufacturing company with strong brand recognition across the globe.

Sylab Scientific delivers high-quality solutions by leveraging technical expertise and components sourced from Elite Thermal Systems Ltd, UK.

The TMH15(Z3)-SY furnace is a bench mounted tube furnace ideal for most general laboratory thermal processing applications.

The TMH15(Z3)/50-SY Three Zone Tube furnace has an internal tube diameter of 50 mm. TMH15(Z3)/50-SY range offers a choice of two heated zone lengths of 450mm & 610mm.

TMH15(Z3)/50/610-SY comes with a heated zone length of 610 mm. Standard zone sizes: 165 mm (zone 1); 290 mm (zone 2); 165 mm (Zone 3). Other zone sizes can also be supplied.



*Illustrative image only

TMH15(Z3)/50/610-SY

SPECIFICATIONS

Maximum Temperature: 1500°C

Maximum Continuous Temperature: 1450°C

- | The TMH(Z3)-SY furnaces are designed to give a longer uniform centre zone temperature than that of the single zone tube furnace models
- | All TMH(Z3)-SY models are controlled by retransmission of set point from the centre zone controller to the end zone controllers
- | This system provides a longer uniform zone temperature than that achieved by the use of single zone furnace of the same length
- | Protective outer mesh covers for improved operator safety
- | Top grade silicon carbide rods in all 3 heated zones
- | Energy efficient, high quality, low thermal mass insulation
- | A work tube is not supplied as an integral part of the furnace and therefore needs to be ordered with the furnace as it is an essential accessory
- | ELITE PID controller offering 8 programs with 8 segments each
- | Temperature Sensor: "R" type Thermocouple
- | Horizontal models are supplied as standard with controls in the base

External Dimensions (mm):	H x W x D	Net Wt.
(Nominal) TMH15(Z3)/50/610-SY	715 x 1000 x 460	65 kg

| **Power Supply:** 400V– 3 Phase

| **Energy rating:** 5.0 kW (Average)
23.0 kW (Peak)

OPTIONS

- | Work tubes of various materials, lengths and diameters for use in the furnace
- | Two or more zones are also available
- | Independent control of each zone is also available
- | The work tubes are available for containment of atmosphere or protection against process contaminants
- | A variety of triple flange gas tight end seals for work tubes to allow processing under Vacuum or gas atmospheres
- | Over Temperature protection controller
- | Multi segment, multi program storage controllers

Vertical versions are available at extra cost, which will have separate temperature console

SYLAB SCIENTIFIC PVT. LTD.

Ground Floor, B-50, Industrial Estate, Sanath Nagar, Hyderabad-500018. INDIA,
Phone: +91 40 67216376 E-mail: sylab@sylabscientific.com Website: www.sylabscientific.com