

# Data Sheet

## Model: TSH15/25/610-SY

1500°C Maximum Temperature,  
Single Zone Horizontal Tube Furnace



### INTRODUCTION

Sylab Scientific Pvt. Ltd., under technical collaboration with Elite Thermal Systems Ltd, UK, manufactures high quality furnaces in India.

Elite Thermal Systems, based in UK, is a 30-year-old furnace manufacturing company with strong brand recognition across the globe.

Sylab Scientific delivers high-quality solutions by leveraging technical expertise and components sourced from Elite Thermal Systems Ltd, UK.

The TSH15-SY furnace is a bench mounted tube furnace ideal for most general laboratory thermal processing applications.

The TSH15/25-SY Single Zone Tube furnace can accept a work tube with internal tube diameter of 25 mm. TSH15/25-SY range offers a choice of three heated zone lengths of 180 mm, 450 mm & 610 mm.

TSH15/25/610-SY comes with a heated zone length of 610 mm.



**TSH15/25/610-SY**

### SPECIFICATIONS

**Maximum Temperature:** 1500°C

**Maximum Continuous Temperature:** 1450°C

| These models are heated by axially mounted silicon carbide elements around the work tube for fast heat up and excellent temperature uniformity

| This furnace design requires the use of a separate work tube of a grade suitable for the maximum temperature rating of 1500°C

| A work tube is not supplied as an integral part of the furnace and therefore needs to be ordered with the furnace as it is an essential accessory

| Temperature Sensor: 'R' type Thermocouple

| Energy efficient, high quality, low thermal mass insulation

| ELITE PID controller offering 8 programs with 8 segments each

| Horizontal models are supplied as standard with controls in the base

|                                  |                  |                |
|----------------------------------|------------------|----------------|
| <b>External Dimensions (mm):</b> | <b>H x W x D</b> | <b>Net Wt.</b> |
| (Nominal) TSH15/25/610-SY        | 715 x 1150 x 460 | 75 kg          |

| Supply: 230V– 1 Phase – 50/60Hz

| Energy rating: 6.0 kW (Average)  
10.1 kW (Peak)

### OPTIONS

| Work tubes of various materials, lengths and diameters for use in the furnace

| The work tubes are available for containment of atmosphere or protection against process contaminants

| A variety of triple flange gas tight end seals for work tubes to allow processing under vacuum or gas atmospheres

| Over Temperature protection controller

| Multi segment, multi program storage controllers

| Thermal Radiation Plugs

**Vertical versions are available which have a separate console**

**SYLAB SCIENTIFIC PVT. LTD.**

Formerly SYLAB PVT. LTD.

Ground Floor, B-50, Industrial Estate, Sanath Nagar, Hyderabad-500018. INDIA,  
Phone: +91 40 67216376 E-mail: sales@sylabscientific.com Website: www.sylabscientific.com