

# Data Sheet

## Model: TSH16/25/610-SY

1600°C Maximum Temperature,  
Single Zone Horizontal Tube Furnace



### INTRODUCTION

Sylab Scientific Pvt. Ltd., under technical collaboration with Elite Thermal Systems Ltd, UK, manufactures high quality furnaces in India.

Elite Thermal Systems, based in UK, is a 30-year-old furnace manufacturing company with strong brand recognition across the globe.

Sylab Scientific delivers high-quality solutions by leveraging technical expertise and components sourced from Elite Thermal Systems Ltd, UK.

The TSH16-SY furnace is a bench mounted tube furnace ideal for most general laboratory thermal processing applications.

The TSH16/25-SY Single Zone Tube furnace can accept a work tube with internal tube diameter of 25 mm. TSH16/25-SY range offers a choice of three heated zone lengths of 180 mm, 450 mm & 610 mm.

TSH16/25/610-SY comes with a heated zone length of 610 mm.



### TSH16/25/610-SY

**External Dimensions (mm):** H x W x D    **Net Wt.**  
(Nominal) TSH16/25/610-SY 825 x 1150 x 562    75 kg

**Power Supply:** 400V– 3 Phase – 50/60Hz

**Energy rating:** 7.0 kW (Average)  
29.0 kW (Peak)

### SPECIFICATIONS

**Maximum Temperature:** 1600°C

**Maximum Continuous Temperature:** 1550°C

- | These models are heated by axially mounted silicon carbide elements around the work tube for fast heat up and excellent temperature uniformity
- | This furnace design requires the use of a separate work tube of a grade suitable for the maximum temperature rating of 1600°C
- | A work tube is not supplied as an integral part of the furnace and therefore needs to be ordered with the furnace as it is an essential accessory
- | Temperature Sensor: 'R' type Thermocouple
- | Energy efficient, high quality, low thermal mass insulation
- | High end Microprocessor PID controller
- | Horizontal models are supplied as standard with controls in the base

### OPTIONS

- | Work tubes of various materials, lengths and diameters for use in the furnace
- | The work tubes are available for containment of atmosphere or protection against process contaminants
- | A variety of triple flange gas tight end seals for work tubes to allow processing under vacuum or gas atmospheres
- | Multi segment, multi program storage controllers
- | Over temperature protection controller
- | Thermal Radiation Plugs

**Vertical versions are available which have a separate console**

### SYLAB SCIENTIFIC PVT. LTD.

Ground Floor, B-50, Industrial Estate, Sanath Nagar, Hyderabad-500018. INDIA,  
Phone: +91 40 67216376 E-mail: sylab@sylabscientific.com Website: www.sylabscientific.com