

# Data Sheet

## Model: TSV12/100/940-SY

1200°C Maximum Temperature,  
Single Zone Vertical Tube Furnace



### INTRODUCTION

Sylab Scientific Pvt. Ltd., under technical collaboration with Elite Thermal Systems Ltd, UK, manufactures high quality furnaces in India.

Elite Thermal Systems, based in UK, is a 30-year-old furnace manufacturing company with strong brand recognition across the globe.

Sylab Scientific delivers high-quality solutions by leveraging technical expertise and components sourced from Elite Thermal Systems Ltd, UK.

The TSV12-SY furnace is a vertical-mounted tube furnace with a separate control console, ideal for most general laboratory thermal processing applications.

The TSV12/100-SY Single Zone Tube furnace has an internal tube diameter of 100mm.

TSV12/100/940-SY comes with a heated zone length of 940 mm.



### SPECIFICATIONS

**Maximum Temperature:** 1200°C

**Maximum Continuous Temperature:** 1150°C

- | The furnace design incorporates an integral elemental tube
- | Protective outer mesh covers for improved operator safety
- | High grade resistance wire wound heating elements
- | Energy efficient, high quality, low thermal mass insulation
- | For aggressive processes, a separate work tube is recommended to minimise the risk of contaminating the elemental work tube
- | A rugged metal sheathed control thermocouple is protected from accidental damage and allows full use of work tube bore
- | Temperature Sensor: 'N' Type thermocouple
- | ELITE PID controller offering 8 programs with 8 segments each
- | **Power Supply:** 230V – 1 Phase – 50/60Hz
- | **Energy rating:** 4.3 kW.

## TSV12/100/940-SY

### OPTIONS

- | Work tubes of various materials, lengths and diameters for use in the furnace
- | The work tubes are available for containment of atmosphere or protection against process contaminants
- | A variety of triple flange gas tight end seals for work tubes to allow processing under vacuum or gas atmospheres
- | Multi segment, multi program storage controller
- | Over temperature protection controller
- | A wide range of furnace mounting configurations, such as Tripod stands, Wall-Mounted Brackets, L-stands, Custom Adjustable Stands, etc.

### SYLAB SCIENTIFIC PVT. LTD.

Ground Floor, B-50, Industrial Estate, Sanath Nagar, Hyderabad-500018. INDIA,  
Phone: +91 40 67216376 E-mail: [sylab@sylabscientific.com](mailto:sylab@sylabscientific.com) Website: [www.sylabscientific.com](http://www.sylabscientific.com)