

Data Sheet

Model: TSV17/75/450

1700 °C Maximum Temperature,
Single Zone Vertical Tube Furnace



INTRODUCTION

The TSV17 furnace is a vertical tube furnace with a separate control panel for most common laboratory thermal processing applications.

The TSV17/75 Single Zone Tube furnace can accept a work tube with internal tube diameter of 75 mm. TSH17/75 range offers a choice of three heated zone lengths of 300 mm, 450 mm and 600mm.

TSV17/75/450 comes with a heated zone length of 450 mm.

SPECIFICATIONS

Maximum Temperature: 1700 °C

Maximum Continuous Temperature: 1650 °C

| These models are heated by substantial molybdenum disilicide loop elements on two sides of the heated chamber.

| This furnace design requires the use of a separate work tube of a grade suitable for the maximum temperature rating of 1700 °C.

| A work tube is not supplied as an integral part of the furnace and therefore needs to be ordered with the furnace as it is an essential accessory.

| Temperature sensor: 'B' type Thermocouple.

| Energy efficient, high quality, low thermal mass insulation.

| High end Microprocessor PID controller.

| Over temperature protection controller.

| Vertical models are supplied as separate control console remote from furnace, separated by 2 meters of flexible interconnecting cable.

| Vertical mounting stand for the furnace body.

| Clamps to hold work tube.



TSV17/75/450

*The above image is for representation only

| Supply: 230V– 1 Phase – 50/60Hz

| Energy rating: 5.5 kW

PHASE CONTROL

Phase angle fired thyristor units operating in conjunction with low voltage secondary isolating transformers providing the correct parameters for the heating elements.

OPTIONAL FEATURES

| Work tubes of various materials, lengths and diameters for use in the furnace.

| The work tubes are available for containment of atmosphere or protection against process contaminants.

| A variety of triple flange gas tight end seals for work tubes to allow processing under gas atmospheres.

| Multi segment, multi program storage controllers.

| Thermal Radiation Plugs.

SYLAB SCIENTIFIC PVT. LTD.

Formerly SYLAB PVT. LTD.

Ground Floor, B-50, Industrial Estate, Sanath Nagar, Hyderabad-500018. INDIA,
Phone: +91 40 67216376 E-mail: sales@sylabscientific.com Website: www.sylabscientific.com