

Data Sheet

Model: TSV18/75/300

1800 °C Maximum Temperature,
Single Zone Vertical Tube Furnace



INTRODUCTION

The TSV18 furnace is a vertical tube furnace with a separate control panel for most common laboratory thermal processing applications.

The TSV18/75 Single Zone Tube furnace can accept a work tube with internal tube diameter of 75 mm. TSV18/75 range offers a choice of three heated zone lengths of 300mm, 450mm & 600mm.

TSV18/75/300 comes with a heated zone length of 300mm.



TSV18/75/300

SPECIFICATIONS

Maximum Temperature: 1800 °C

Maximum Continuous Temperature: 1750 °C

- | These models are heated on both sides of the chamber by molybdenum Tungsten Disilicide elements.
- | This furnace design requires the use of a separate work tube of a grade suitable for the maximum temperature rating of 1800 °C.
- | A work tube is not supplied as an integral part of the furnace and therefore needs to be ordered with the furnace as it is an essential accessory.
- | Temperature Sensor: 'B' Type Thermocouple.
- | Energy efficient, high quality, low thermal mass insulation.
- | High end Microprocessor PID controller.
- | Over temperature protection controller.
- | Vertical models are supplied as separate control console remote from furnace, separated by 2 meters of flexible interconnecting cable.
- | Vertical mounting stand for the furnace body.
- | Clamps to hold work tube.

*The above image is for representation only

- | Supply: 230V– 1 Phase – 50/60Hz
- | Energy rating: 5.0 kW

PHASE CONTROL

Phase angle fired thyristor units operating in conjunction with low voltage secondary isolating transformers providing the correct parameters for the heating elements.

OPTIONAL FEATURES

- | Work tubes of various materials, lengths and diameters for use in the furnace.
- | The work tubes are available for containment of atmosphere or protection against process contaminants.
- | A variety of triple flange gas tight end seals for work tubes to allow processing under gas atmospheres.
- | Multi segment, multi program storage controllers.
- | Thermal Radiation Plugs.

SYLAB SCIENTIFIC PVT. LTD.

Formerly SYLAB PVT. LTD.

Ground Floor, B-50, Industrial Estate, Sanath Nagar, Hyderabad-500018. INDIA,
Phone: +91 40 67216376 E-mail: sales@sylabscientific.com Website: www.sylabscientific.com